# AMPLATZER® Exchange System

## Instructions for Use

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STERILE EO

600071-003

Amplatzer is a registered trademark of AGA Medical Corporation  
Not in any way connected with medical gas or equipment sold under the "AGA" brand by AGA AB or its successors.

## AMPLATZER® EXCHANGE SYSTEM Instructions For Use

### Intended Use


The AMPLATZER® Exchange System is intended for removal of an AMPLAZER Delivery Sheath and subsequent exchange for an Amplatzer Delivery Sheath of equal or larger diameter.

### Description

The AMPLATZER Exchange System or "Bail out" mechanism is a delivery system especially adapted for use in conjunction with the AMPLATZER family of occlusion devices. The system components are identical to the AMPLATZER Delivery System, with the exception of the dilator, which incorporates an enlarged inner lumen for passage over an AMPLATZER delivery cable.

The AMPLATZER Exchange System incorporates the following elements:

- Delivery sheath with Touhy-Borst adapter
- Dilator – used to ease penetration of tissue.
- Loading Device – used to introduce an AMPLATZER occluder into the delivery sheath.
- Plastic Vise – facilitates direction control and serves as the "handle" for disconnecting (unscrewing) the delivery cable from the device.
- Delivery Cable – the device is screwed onto the distal tip of the delivery cable, which allows for placement (and if necessary, retrieval) of the device.

- Delivery Sheath** - The delivery sheath is manufactured from PTFE\*, a material that is widely used in many medical devices. The delivery sheath allows for delivery of the AMPLATZER device. A hemostasis valve prevents blood reflux and air aspiration
- Dilator** - A locking tapered dilator is added for ease of tissue penetration and is a removable component. The dilator is manufactured from PTFE and has a lumen large enough to accommodate the thickness of the AMPLATZER delivery cable.
- Loading Device** - Also known as the Introducer, this component is used to introduce the AMPLATZER Occluder into the delivery sheath. The Introducer is also manufactured from PTFE.
- Delivery Cable** - The shaft of the delivery cable is constructed from stainless steel in a twisted configuration. An outer coil prevents kinking and also increases push strength. The outer surface is uncoated. A "mini" threaded screw is welded to the distal end for attachment to the AMPLATZER device. A thread on capsule located on the proximal end of the cable permits connection of two delivery cables thus doubling the length of the system.
- Plastic Vise** - Also referred to as a pin vise, this standard accessory item is manufactured from polycarbonate material. This is a non-tissue contact component. The pin vise is utilized to facilitate directional control and serves as the "handle" for disconnection (unscrewing) of the delivery cable from the AMPLATZER device.

### Warnings

Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse reactions if reused. Accordingly, AGA Medical Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.

### Precautions

- Store in a cool, dry place.
- Do not use if package is open or damaged.
- Use prior to "Use By" date on the label.
- Do not expose to organic solvents.
- Do not autoclave. Exposure to temperatures above 54°C (130°F) may damage the components.
- If strong resistance is met during manipulation, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of the resistance cannot be determined, withdraw all components.

\*PTFE is an acronym for the chemical compound polytetrafluoroethylene.

### Complications

Procedures requiring percutaneous catheter/guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Possible complications include, but are not limited to:

- Air embolism
- Hematoma at the puncture site
- Infection
- Perforation of the vessel wall.

### Procedure

- Peel open the package and flush the delivery sheath, dilator and loader using sterile heparinized normal saline or a similar isotonic solution by connecting a syringe to the luer hub of each component.
- Insert the dilator into the delivery sheath and secure with the locking mechanism.
- Screw the delivery cable to the capsule located on the proximal end of the original delivery cable thus doubling its length.

- Once securely attached, remove the original delivery sheath from the patient and discard.
- Introduce the dilator/delivery sheath exchange assembly over the lengthened cable assembly through the groin. Once the delivery sheath has reached the inferior vena cava, remove the dilator to allow back bleeding to purge all air from the system.
- Advance the delivery sheath over the cable through the communication and disconnect the secondary cable.
- Proceed with device recapture and/or deployment according to the instructions provided with the AMPLATZER Occlusion device.

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